



Paediatric catalogue

Only the best is good enough for our children...

Dear Customer,

At DJO we are dedicated to producing the very best hi-tech orthopaedic products, through intensive research, clinical trials, and award-winning scientific design.

For healthcare professionals working in Paediatric departments and those regularly consulting infants and children, our new catalogue sets out to provide and showcase DJO's products that are available for this important patient group.

The catalogue not only illustrates children's products from our DonJoy, Procure and Aircast ranges but also includes Axmed's* well-established paediatric products which have recently been integrated into DJO's portfolio under the DonJoy brandname.

The collection of products presented in this catalogue combine those specifically designed for children and only available in children's sizes, and the smaller sizes from our key adult product lines. All products are CE marked in line with current EU regulations.

We hope this catalogue will help increase awareness of the products that exist to provide our children with the optimal care they deserve.

DJO is a company that aims to "never stop getting better" so we are always interested in receiving suggestions to help us develop our products, to help improve your patients' comfort and accelerate their recovery. We always consider suggestions, so please either complete and return the enclosed customer feedback reply form or email me personally at lieve.vandenbergh@djortho.com

Finally, please do visit our new International website, www.djo.eu, for our latest news, research updates and product information and also the opportunity to subscribe to our new quarterly newsletter.

Yours sincerely,



Dr. Lieve Vanden Berghe, PT PhD
DJO, Director of International Marketing



*Axmed is a French soft goods company that was acquired by DJO in 2006. The Axmed brand is abandoned and the quality products are now sold under the DonJoy brand name

Cervical Collars



Part Number	79-83130
Height Adj.	1.5", 1.75", 2" (4cm, 4.5cm, 5cm)
Neck Circ.	8"-18" (20.5cm-46cm)
Size	Paediatric

Patriot® Collar

The Patriot® Collar is a one-piece injection molded collar which provides proper support and helps restrict cervical spine flexion, extension and rotation to promote patient recovery following trauma. Constructed with durable high-density polyethylene, lined with a soft closed cell foam for maximum patient comfort. The collar has multiple heights to fit most adults and paediatrics. An oversized trachea opening provides quick access to carotid pulse monitoring and emergency tracheotomies. Stacks flat to preserve space and is X-Ray, CT and MRI lucent.

Does not contain Latex.



Part Number	79-83161
Height	2.25" (6cm)
Size	Paed, XS

Height = Sternal notch to chin.

California™ Cervical Collar

The California™ Cervical Collar is a pre-formed, lightweight, two-piece polyethylene foam support, including rigid plastic supports at the front and rear. The collar maintains the head in a neutral position and may be used to provide support in combination with a cervical or scoliosis cast. Ideal for cervical spine immobilisation and extraction.

Does not contain Latex.



Part Number	79-83011
Height x Length	2" x 14" (5cm x 36cm)
Neck Circ.	7" x 12" (18cm x 30.5cm)
Size	Paed/XXS

Form Fit Cervical Collar

The Form Fit Cervical Collar is a medium density foam collar which helps support the cervical spine in a neutral position. It has a contour design for use in either a flexed or extended position. With a cotton stockinette cover for better patient hygiene, the Form Fit is ideal for post-operative rehabilitation and cervical sprains and strains.

Does not contain Latex.



Contoured Clavicle

The Contoured Clavicle has pre-formed straps for maximum comfort. It is constructed of bonded poly-urethane foam and covered with a cotton stockinette. Featuring a figure 8 design with a 2-way contact closure or 3-way buckle closure available. Ideal for clavicular fractures and postural problems.



Does not contain Latex.

Part Number 79-98671
Size Paediatric

Part Number 79-98672
Size XXS

Clavicle Splint

The Clavicle Splint has adjustable prong buckles which ensure a correct fit and prevents slippage. With stockinette covered felt shoulder straps, the Clavicle Splint is washable and ideal for clavicular fractures and postural problems.



Does not contain Latex.

Part Number 79-85001
Circ. 18"-20"
(46cm-51cm)
Size Paed, XXS

Clavicle Supports



Arm Slings & Shoulder Immobilisers



Part Number	79-99132
Depth x Length	4.5" x 11" (11.5cm x 28cm)
Size	Paed/XS
Part Number	79-99133
Depth x Length	6.5" x 13" (11.5cm x 28cm)
Size	Child/S

Part Number	79-84001
Depth x Length	6.5" x 8.25" (16.5cm x 21cm)
Size	Paed/XXS
Part Number	79-84002
Depth x Length	7" x 11" (18cm x 28cm)
Size	Child/XS

Part Number	226-P
Length of Forearm	<9" (<23cm)
Size	Paediatric

Healthcare Bear Arm Sling

This arm sling has a cotton/poly envelope with paediatric print and a cotton web strap with contact closure for easy application. Designed to help ensure immobilisation of the shoulder after injury in infants. Can be used in conjunction with a cast.



Does not contain Latex.

Deluxe Arm Sling W/Pad

The Deluxe Arm Sling has a generous size envelope to support the arm and a cotton/poly construction with web strap and contact closure. A thumb loop helps prevent the sling from migrating. Shoulder pad included. Ideal for cast support; injuries to the arm, wrist or hand.



Does not contain Latex.

IMMO

The Immo Sling & Swathe is recommended for immobilisation after upper extremity trauma such as post-surgical repair and humeral lesions. Also indicated for paralysis and hemiplegia.



Does not contain Latex.

Cock-Up Splint



The Cock-up Splint is a durable, canvas construction with malleable aluminum palmer stay. Lace closure for ease of application. Ideal for post-fracture casting, wrist sprains/strains, Carpel Tunnel Syndrome.
Does not contain Latex.

Right	Left	Length	Size
79-87321	79-87331	4.5" (11.5cm)	Paed/XXXS
79-87341	79-87351	6" (15cm)	Paed/XXS
79-87342	79-87352	6" (15cm)	Paed/XS
79-87381	79-87391	8" (20.5cm)	Paed/XXS
79-87382	79-87392	8" (20.5cm)	Paed/XS

Wrist circumference

XXXS= up to 3.5" (9cm)
XXS=3.5"-4.5" (9cm-11.5cm)
XS=4.5"-5.5" (11.5cm-14cm)

Foam Wrist Splint



This wrist splint has a foam tricot construction for breathable comfort. Dorsal and palmer malleable aluminum stays may be removed and contoured to proper angulation. Ideal for Carpel Tunnel Syndrome and immobilisation of the wrist.

Does not contain Latex.

Right	Left
Part Number 79-87422	Part Number 79-87432
Wrist Circ. 4.5" - 5.5" (11.5cm-14cm)	Wrist Circ. 4.5" - 5.5" (11.5cm-14cm)
Size Paed/XS	Size Paed/XS

Loop-Lock Cockup Splint



The Loop-Lock Cockup Splint has a malleable aluminum palmer stay for custom anatomical position and a sturdy canvas construction with loop-lock closure. The splint allows full finger function above the palmer crease and is ideal for Carpel Tunnel Syndrome and immobilisation of the wrist.
Does not contain Latex.

Right	Left
Part Number 6" (15cm) 79-87362	Part Number 6" (15cm) 79-87372
Part Number 8" (20.5cm) 79-87402	Part Number 8" (20.5cm) 79-87412
Size Paed/XS	Size Paed/XS

Wrist, Forearm & Elbow Supports



Wrist, Forearm & Elbow Supports



Universal Paediatric Wrist Splint

Constructed of a breathable foam laminate, the Universal Paediatric Wrist has a stockinette added at the strap area for additional patient comfort. Sizes include paediatric and child. Indications include: Sprains and strains, Carpal Tunnel Syndrome and post-operative rehabilitation for the wrist.



Does not contain Latex.

Right		
Part Number	Length	Size
79-87420	4.5" (11.5cm)	Paed
79-87440	6" (15cm)	Child
Left		
Part Number	Length	Size
79-87430	4.5" (11.5cm)	Paed
79-87450	6" (15cm)	Child

Elbow Immobiliser

The Elbow Immobiliser has a soft foam/nylon laminate construction for patient comfort. Aluminum stays help prevent patients from bending the elbow joint or accessing tubes. The immobiliser also has contact closure straps for ease of application. Indications include: limb positioning, post-operative cleft palate and IV therapy. Does not contain Latex.



Part Number	Length	Width	Size
79-91512	4" (10cm)	7.5" (11cm)	Infant/XS
79-91513	6" (15cm)	9.5" (24cm)	Child/S
79-91515	7" (18cm)	12" (30.5cm)	Youth/Adult/M

Colles Splint

The Colles Splint is a premolded aluminum splint designed for wrist support and immobilisation prior to surgery or casting. Ventilation holes provide air circulation for patient comfort.



Does not contain Latex.

Padded		Length x Width	Size
Right	Left		
79-72122	79-72112	5.75" x 2.5" (14.5cm x 6.5cm)	Paed/XS
Unpadded		Length x Width	Size
Right	Left		
79-72142	79-72132	5.75" x 2.5" (14.5cm x 6.5cm)	Paed/XS



AT4

The AT4 is a one-piece foam construction with breathable Waylex™ mesh lining and includes 2 removable double bend posterior stays for immobilisation of the knee in slight flexion. It has a buckle closure for easy application and fits either the left or right knee. Ideal for post-operative total immobilisation after ligament injury or other knee joint injury.

Part Number	AT430/01
Length	11.8" (30cm)
Size	Paediatric

Does not contain Latex.



Super Knee Splint

The Super Knee Splint is a comfortable, sturdy, foam/mesh construction with alternating web straps and loop/lock contact closure. Medial/lateral I-beam stays and three double bend posterior stays provide proper immobilisation. Adjustable foam pop pad included. Ideal for immobilisation of the knee joint.

Part Number	79-80002
Circ.	12"-14" (30.5cm-35.5cm)
Length	12" (30.5cm)
Size	Paed/XS

Does not contain Latex.

Knee Immobilisers



Knee Supports



Each Defiance III® is custom made


	Minimum Size
Thigh:	12.5" (32cm)
Calf:	10.25" (26cm)

Defiance III®

The new Defiance III® Knee Brace has taken the best from our Defiance® brace platform and added a "host of options" for our customers. The enhancements are ideal for skiing and riding type activities (Motocross and Mountain Biking), bilateral bracing and for difficult to fit leg shapes. This brace provides superior ACL, PCL, MCL and LCL instability support for all activity levels.

- Patented Four-Points-of-Leverage™ design
- ACL, PCL and CI versions
- Swooping medial thigh frame for increased inner thigh clearance
- Swivelling strap tabs to match anatomy through movement
- Anti-Migration calf strap and Supra Condylar Pad for brace suspension
- FourcePoint™ and Standard polycentric hinges available
- Lightweight, low profile design
- Short calf version available
- Assortment of colour, graphic and gradation options available
- Lifetime limited warranty on frame and hinge

Does not contain Latex.



ACL Everyday

The ACL Everyday represents DonJoy's most comfortable, easiest to use, patient-ready knee brace and is ideally suited for non-contact sports and activities of daily living. Recommended for mild ACL/MCL/LCL instabilities, offering comfortable support for the activities of daily living or non-contact sports.

- Total Comfort Guarantee
- Four-Points-of-Leverage™ system pioneered by DonJoy®
- Lightweight (15 oz.)
- Sleek, low profile (.080" thick), aesthetically pleasing, aircraft grade aluminum frame

Does not contain Latex.



Part Number	11-1601-1
Leg	Right
Size	XS / 13"-15.5" (33cm-39.5cm)

Part Number	11-1600-1
Leg	Left
Size	XS / 13"-15.5" (33cm-39.5cm)

Length 14" (36cm)

Measurements taken
15cm above mid-patella



Air-Stirrup® Ankle Brace

The Air-Stirrup® has anatomically designed semi-rigid shells lined with the patented Duplex™ aircell system for support and graduated compression during ambulation. The compression promotes efficient edema reduction in addition to helping accelerate rehabilitation from ankle sprains.

Does not contain Latex.

Right
Part Number 02JR
Height 6" (15cm)
Size Paediatric

Left
Part Number 02JL
Height 6" (15cm)
Size Paediatric



Surround® FLOAM™ Ankle

This brace provides stability to the ankle by minimising ankle rotation. Comprised of two polypropylene outer shells connected with a soft, adjustable foot piece. A unique FLOAM™ bladder provides uniform compression and a customised fit. FLOAM™ has no memory and therefore prohibits pressure from being displaced to bony prominences. Fits right or left foot. Ideal for acute ankle sprains, chronic ankle instabilities and osteoarthritis.

Does not contain Latex.

Part Number 81-81193
Height 6" (15cm)
Size Paediatric

Ankle Braces




Post-Op and Cast Shoes



Silicone Heel Cups

Part Number 79-81102
Shoe Size XS

The Silicone Heel Cups are molded medical-grade soft silicone to help reduce shock and vibration loads to joints, to address plantar fasciitis, heel spurs and Achilles tendonitis. They have a softer center in the heel region to help reduce strain to the area sensitive to pain. Sold in pairs.
Does not contain Latex.



Part Number 79-81110
Length 4.5" (11.5cm)
Toe Width 2.5" (6.5cm)
Size XXXS

Part Number 79-81111
Length 5.5" (14cm)
Toe Width 3.0" (7.5cm)
Size XXS

Rocker Cast Boot

The Rocker Cast Boot has a rubber rocker sole to allow for more natural gait during ambulation. It is a sandal style boot with a durable canvas upper and forefoot with heel loop-lock contact closure for ease of application. For use with lower leg casted fractures.

Does not contain Latex.



Part Number 79-81222
Length 7.25" (18.5cm)
Toe Width 3.5" (9cm)
Size Paediatric

Rocker Cast Shoe w/Flex Sole

Constructed of a durable canvas upper with a flexible rocker sole for more natural gait during ambulation, the Rocker Cast Shoe has a contact closure forefoot and heel straps for ease of application. Ideal for lower leg casted fractures.

Does not contain Latex.



Part Number 79-90192
Childs Shoe Size 13.5-2.5 (32.2-35)
Size XS

Part Number 79-90193
Adult Shoe Size 3.5-5.5 (35.5-38)
Size S

Post-Op Shoe

The Post-Op Show has a rigid E.V.A. rocker sole and contoured heel for added patient comfort. It includes a low-profile foam/nylon mesh upper with tricot lining and a protective tongue provides maximum comfort. Ideal for post foot surgery, metatarsal injuries and soft tissue procedures.

Does not contain Latex.



Paediatric Male
Part Number 79-81132
Shoe Size 4.5-6.5 (37-39)
Size Paed/XS

Paediatric Female
Part Number 79-81142
Shoe Size 13.5-2.5 (32.2-35)
Size Paed/XS

Medical/Surgical Shoe

Designed for a variety of post-op and post-trauma applications, this shoe has a soft foam mesh upper with padded tricot lining. It features a semi-flexible skid resistant sole and comfortable contoured heel with a large overlapping toe area to accommodate for post-operative dressing and swelling. Fits right or left foot.

Does not contain Latex.



Aircast® XP Walker™



The XP Walker™, the next generation of the Aircast® pneumatic walking brace, has new features to increase comfort, simplify use, and promote compliance while still retaining the clinically effective elements found in the original prefabricated walking brace by Aircast®. Each XP Walker™ features a low rocker sole for added comfort and ease of ambulation, a wider foot base ensuring plenty of room for dressings without sacrificing comfort and relocated valves to improve ease of use and patient compliance. The XP Walker™ can be used in conjunction with Aircast® Heel Supports as Achilles Wedges.

Does not contain Latex.

Part Number 01P-P
Adult Shoe Size <4 (<36)

Aircast® FP Walker™



The next generation of the cost-effective Aircast® foam pneumatic walking brace, the FP Walker™ provides full-shell protection with two distal aircells for pneumatic support. Each FP Walker™ has a low rocker sole for added comfort and ease of ambulation and a wider foot base with plenty of room for dressings without sacrificing comfort.

Does not contain Latex.

Part Number 01F-P
Adult Shoe Size <4 (<36)



Walkers

Walkers



Part Number 11-1382-1-00000
Shoe Size <2 (<34)

Part Number 11-1382-2-00000
Shoe Size 2-3 (34-35)

Maxtrax® Walker

The MaxTrax® is designed with a wider foot bed for comfort and stability. The low profile, rocker bottom helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. The walker also features a washable, one-piece liner which has been designed to keep seams away from post-surgical areas. The Maxtrax® walker is recommended for use following stress fractures of the lower leg, soft tissue injuries (grade 2 and 3 sprains), stable fractures and injuries of the foot and ankle, bunionectomies and metatarsal fractures.

Does not contain Latex.



Part Number 11-1371-1
Shoe Size <2 (<34)

Maxtrax® Air Walker

The MaxTrax® Air walkers are designed to increase stability and help decrease pain and edema following trauma or post-operative procedures. They have cushioned inner/outer soles designed to absorb shock upon heel strike and improve patient comfort. An integrated pump makes the liner simple to inflate and adjust and a low profile, rocker bottom helps promote natural gait and reduce plantar pressures.

Does not contain Latex.



Paediatric Ankle Cryo/Cuff™



The Paediatric Ankle Cryo/Cuff™ is specifically designed to treat children's injuries (ages 1 to 7). Using a combination of controlled compression and cold, the Cryo/Cuff™ minimises haemarthrosis, swelling and pain.

Does not contain Latex.

Part Number | 10P01

Paediatric Knee/Elbow Cryo/Cuff™

The Paediatric Knee/Elbow Cryo/Cuff™ is specifically designed to treat children's injuries (ages 1 to 7). It combines the therapeutic benefits of controlled cold compression to minimise haemarthrosis and swelling and reduce pain. Due to the unique design of the cuff, it can be utilised effectively on both the knee and elbow of the patient.



Does not contain Latex.

Part Number | 11P01

Cryo Therapy



Child Care, Safety & Traction

Part Number	79-90512
Chest Circ.	11.5" - 14.5" (29cm - 37cm)
Size	Paed, XS

Pavlik Harness

Specifically designed for treatment of hip dysplasia and dislocation, the Pavlik Harness has a chest band and shoulder straps made with heavy webbing with Kodel Polyester pile backing for durability and comfort. Combining felt and web stirrups with contact closure for ease of application.

Does not contain Latex.



Part Number	79-90113
	1 pair
Size	Paed, XS

Foam Limb Restraints

The Foam Limb Restraints have contact closure with double "D" ring cinch lock. Available in pairs. Indications include: positioning of limbs and limb restraints.

Does not contain Latex.



Part Number	79-90122
	10 pairs
Size	Paed, XS

Double Buckle Limb Restraint

The Double Buckle Limb is a foam restraint for the ankle or elbow with double-buckle closure and an extra long tie. Sold 10 pairs per pack. Indications include: positioning of limbs and limb restraints.

Does not contain Latex.



Asia-Pacific:

DJO Asia-Pacific Corporate Office
Aircast Asia-Pacific Limited
Rm. 701, 7/F Tower 1
Grand Central Plaza,
Shatin, Hong Kong
Tel.: +852/ 310 514 15
Fax: +852 / 310 514 44
Email: phoyer@aircstasia.com

Benelux:

DJO Benelux BVBA
NCI Business Center
Twin Squares - Building Vendôme
Culliganlaan 1B
B-1831 Diegem
BELGIUM
Tel.: +32 2 403 1320/19
Fax: +32 2 403 1212
Email: Sarah.vanwijnsberghe@djortho.com

Canada:

DJO Canada
2835 Argentia Road, Unit # 5
Mississauga
ON L5N 8G6
CANADA
Tel.: + 1 1866 866 5031
Fax: + 1 1866 866 5032
Email: Miles.white@djortho.com

France:

DJO France S.A.S.
Z.I. Les Pontôts
8 Allée du Canal
64600 Anglet
FRANCE
Tel.: +33 559 52 86 90
Fax: +33 559 52 86 91
Email: Sce.Cial@djortho.com

Germany:

DJO Deutschland GmbH
Kulmbacher Strasse 51
D-95512 Neudrossenfeld
GERMANY
Tel.: +49 (0)9203 97350
Fax: +49 (0)9203 973519
Email: de.info@djortho.com

Italy:

DJO Italia Srl
Via Leonardo Da Vinci 97
Trezzano Sul Naviglio
20090 Milano
ITALY
Tel.: +39 02 484 63386
Fax: +39 02 484 09217
Email: itinfo@djortho.com

Denmark, Finland,**Norway & Sweden:**

DJO Nordic A/S
Hovedgaden 51, 1. sal
2970 Hørsholm
DENMARK
Tel.: +45 4541 0670
Fax: +45 4541 0671
Email: info.nordic@djortho.com

Spain & Portugal:

DJO Ibérica
Calle Gabiria, 15
20305 Irun
SPAIN
Tel.: +34 943 638 167
Fax: +34 943 638 174
Email: svc.cial@djortho.com

UK & Ireland:

DJO UK Ltd.
7 The Pines Business Park
Broad Street
Guildford
Surrey
GU3 3BH
ENGLAND
Tel.: +44 (0)1483 459 659
Fax: +44 (0)1483 459 470
Email: ukorders@djortho.com

